

Charm Fast Phage Test

for Fecal Indicator in Drinking Water

8
hours

Product Overview



Positive



Negative
<1 pfu/100 mL

Charm Sciences' EPA Approved Fast Phage test detects coliphage, a viral indicator of fecal contamination in water, in less than 8 hours. With its fast results, Charm® Fast Phage enables municipalities and public health officials to be proactive – providing rapid remediation of vital risks in water supplies and wastewater management.

Current tests for coliphage take 36 - 60 hours for result. Fast Phage takes only 16 hours to confirm results and has an 8-hour predictor for fecal contamination, enabling a same-day public health response.

Complement your bacterial testing. Since 1990, viral contamination has been the source of more than half of all of waterborne illnesses^A. Bacterial (coliform/*E. coli*) indicator tests for drinking water are blind to viral contamination, so even if water passes bacterial testing, it could still contain viral risk.

^AOccurrence and Monitoring Document for the Final Ground Water Rule, October 2006
http://www.epa.gov/ogwdw/disinfection/gwr/pdfs/support_gwr_occurrence-monitoring.pdf

Charm Fast Phage Test:

- ✓ **Detects positives in less than 8 hours**
- ✓ **Confirms negatives in 16 hours**
- ✓ **Augments bacterial testing with viral indicator**
- ✓ **Contains pre-measured reagents, unit-dosed, and ready-to-use *E. coli* host**
- ✓ **Delivers clear, easy-to-read fluorescent results**

Simple Procedure^B



Start




**Transfer For Fluorescent
Predictor**
5 hrs



**Transfer For Plaque
Confirmation Assay**
5 hr



Fluorescent Warning
8 hrs




Plaque Confirmation
16 hrs

^B For detailed explanation on procedure, please see Operator's Manual

Positive Prediction

Positive coliphage prediction in a working day enables faster response to fecal risks in the water supply and compliments bacterial test results with a viral indicator.

Response to Fecal Risks in the Water Supply



Ordering Info

Order Codes	Each Kit Includes
FP-SOM-25K FP-FPLUS-25K	25 Phage Somatic Presence/Absence Tests 25 Phage F+ Male Specific Coliphage Presence/Absence Tests

See Operator's Manual for additional kits and order codes