


ROSA FAST Aflatoxin Quantitative Test for Feed and Grain

3 or 5
minutes

Product Overview



The **ROSA® FAST Aflatoxin Quantitative Test (AFQ-FAST)** delivers rapid and accurate results using ROSA (Rapid One Step Assay) lateral flow technology. ROSA mycotoxin tests are simple to run; just add the diluted sample extract to the strip, incubate for 3 or 5 minutes, and read results in the Charm EZ-M system.

Benefits of AFQ-FAST Test Include:

- Range: 0-30 ppb, 20-150 ppb, 50-500 ppb
- Results in 3 or 5 minutes
- FGIS Approval

Run on the Charm EZ-M System:

- High volume output with batching using multiple quad incubators
- Standard curve built into reader reducing testing cost by eliminating end-user calibration
- Printable or downloadable results to any computer or LIMS system

Commodities

FGIS Approved Commodities (70% Methanol)

3 Min	Corn				
5 Min	Barley	Brown Rice	Corn Germ Meal	Corn Gluten Meal	Corn/Soy Blend
	DDG	DDGS	Hominy	Millet	Oats
	Popcorn	Rough Rice	Rye	Sorghum	Soybean Hulls
	Soybean Meal	Soybeans	Wheat		

^AFGIS - Federal Grain Inspection Services

Charm® Validated Commodities (70% Methanol)

5 Min	Milled Rice				
-------	-------------	--	--	--	--

Charm Validated Commodities (84% Acetonitrile)

5 Min	Corn Gluten Feed	Copra Meal	Corn Silage	Cottonseed Hulls	Cottonseed Meal
	Palm Kernel Meal	Peanut Hulls	Peanut Meal	Rice Bran	Whole Cottonseed



Simple Procedure



Add Sample.

Incubate.

Read results in less than 5 seconds!



For batch processing, incubate multiple strips simultaneously, then read results in the Charm EZ-M system.

Ordering Info

Order Codes	Each Kit Includes
LF-AFQ-FAST-20K LF-AFQ-FAST-100K LF-AFQ-FAST-500K	Test Strips, Positive Control, Dilution Buffer

See Operator's Manual for additional kits and order codes