

ROSA® FAST5 Fumonisin Quantitative Test Flow Chart



See Approved Commodities Below

**Test Ranges: 0 to 1.5 ppm
0 to 6 ppm**

Charm-Validated Commodities:

A2:1 Extraction Ratio: Barley, Corn, Corn Gluten Meal, Flaking Corn Grits, Millet, Oats, Rough Rice, Sorghum, Wheat

B3:1 Extraction Ratio: Distillers Dried Grain with Solubles (DDGS)

C3:1 Extraction Ratio: Corn Germ Meal

Sample Preparation



(1) Weigh

Ground sample^D

50 g
or 10 g



(2) Add Solvent


70% Methanol

100 mL^A/ 150 mL^{BC}
or 20 mL^A/30 mL^{BC}



(3) Extract

Shake vigorously or blend for 1 minute; do not exceed 2 minutes



(4) Clarify

Centrifuge, filter or allow sample to settle



(5) Dilute^E


Prepare Diluted Extract

100 µL Extract* + 1.0 mL FUM Dilution Buffer
*158 µL Extract for Corn Germ Meal



Diluted Extract

Filter for:
Barley and Wheat




Pass Diluted Extract through RC15 filter

^D50 g is the Charm recommended sample weight, 10 g is an optional sample weight

^EDiluted Extract for corn germ meal samples cannot be assayed, continue to prepare and assay 2nd Diluted Extract


Test Procedure



(1)

Place test strip in ROSA Incubator or Charm EZ[®]-M system.

For Charm EZ-M system select appropriate test, commodity and dilution if prompted.



(2)

Peel tape.
Pipet 300 µL Diluted Extract into sample compartment.
Reseal tape.




(3)

Close lid.
Incubate for 5 minutes.

Read Result

(1) Inspect test strip




(2) Read result with ROSA-M Reader or Charm EZ-M system

ROSA-M Reader: Select FUM channel in 3-line mode (blinking) and appropriate MATRIX.

Charm EZ-M system: Select appropriate test, commodity and dilution if prompted.

Extraction Ratio	Sample (Dilution)	MATRIX	Quantitation Range	LOD
2:1^A	Diluted Extract (DE)	00	0.5 to 1.5 ppm	0.25 ppm
	2 nd Diluted Extract (2ND DE)	01	1 to 5.4 ppm	-
3:1^B	Diluted Extract (DE)	03	0.5 to 1.5 ppm	0.25 ppm
	2 nd Diluted Extract (2ND DE)	04	1 to 5.4 ppm	-
3:1^C	2 nd Diluted Extract (2ND DE)	01	1 to 5.4 ppm	-

For quantitation of 1 to 5.4 ppm:



(1) Prepare 2nd Diluted Extract

(2) Repeat Test Procedure (steps 1, 2, 3) with 2nd Diluted Extract

(3) Read Result