


Charm QUAD Test for Milk

5
minutes

Product Overview




Charm® Rapid One Step Assay (ROSA®) technology is trusted by the dairy industry worldwide for fast and accurate milk residue screening. The **Charm QUAD® Test (QUAD)** detects chloramphenicol, beta-lactam, tetracycline, and streptomycin drugs in raw commingled milk meeting the Russian Technical Regulation of the Customs Union. The QUAD test may also be used to test at or below EU/CODEX Maximum Residue Limit (MRL) and EU Minimum Required Performance Limit (MRPL).


Benefits of QUAD Test Include:

- Achieves results in 5 minutes
- Detects chloramphenicol, beta-lactam, tetracycline, and streptomycin drugs
- Incubates at 40 ± 1 °C
- Meets Russian Technical Regulation of the Customs Union
- Meets EU/CODEX MRL and EU MRPL
- Validated by ILVO – Institute for Agricultural and Fisheries Research (ILVO), Belgium

Simple Procedure




Add Sample.



Incubate.



Read results in less than 5 seconds!



For batch processing, incubate multiple strips simultaneously, then read results in the Charm EZ system.



Sensitivity in Milk

Sensitivity-Detection in Cow Milk at 0 to 7 °C					
Beta-lactam Drug	Detection Range ^A (ppb ^B)	EU/CODEX MRL (ppb ^B)	Tetracycline Drug	Detection Range ^A (ppb ^B)	Russian Technical Regulation ^D (ppb ^B)
Amoxicillin	3 to 5	4	Chlortetracycline	4 to 6	10
Ampicillin	3 to 4	4	Doxycycline	6 to 10	10
Cefacetile	20 to 40	125	Oxytetracycline	4 to 6	10
Cefalonium	8 to 12	20	Tetracycline	3 to 6	10
Cefazolin	20 to 40	50	Amphenicol Drug	Detection Range ^A (ppb ^B)	EU MRPL ^E and Russian Technical Regulation (ppb ^B)
Cefoperazone	1 to 3	50	Chloramphenicol	0.2 to 0.3	0.3
Cefquinome	10 to 20	20			
Ceftiofur and Metabolites ^C	30 to 40	100 / 100	Aminoglycoside Drug	Detection Range ^A (ppb ^B)	EU/CODEX MRL and Russian Technical Regulation ^D (ppb ^B)
Cephapirin	20 to 35	60	Dihydrostreptomycin	75 to 150	200 / 200
Cloxacillin	20 to 30	30			
Dicloxacillin	20 to 30	30			
Penicillin G	2 to 3	4 / 4	Streptomycin	75 to 150	200 / 200

^A Positive at least 95% of the time.

^B parts per billion, µg/kg or µg/L

^C Ceftiofur parent drug sensitivity is approximately 1/2 that reported in the table.

^D Russian Technical Regulation of the Customs Union on Safety of Milk and Dairy Products (TR TS 033/2013)

^E Minimum Required Performance Limit

Ordering Info

Order Codes	Each Kit Includes
LF-QUAD-RUS-20K LF-QUAD-RUS-100K	Test Strips and Positive Controls

See Operator's Manual for additional kits and order codes



659 Andover Street, Lawrence, MA 01843-1032, USA
T +1.978.687.9200 | **F** +1.978.687.9216 | **E** info@charm.com | **www.charm.com**

© 2019 Charm Sciences, Inc. Charm, Charm EZ, ROSA, and Running Cow image are registered trademarks, and test strip product packaging is trade dress, of Charm Sciences, Inc.

MRK-852-006 Oct-2019