


Charm II : Aflatoxin Analysis of Dairy Products



Charm II 6600/7600 Analyzer

CHARM II AFLATOXIN TEST AT A GLANCE

- ▶▶ Rapid test to meet EC MRLs or US safe levels
- ▶▶ Effectively screens at 25 ppt aflatoxin M1 in milk or 500ppt (.50ppb)
- ▶▶ Multiple samples can be run simultaneously in under 15-30 minutes [depending on assay]
- ▶▶ Uses low cost reagents



CHARM II AFLATOXIN ASSAYS

CII ASSAY	ASSAY TIME	POSITIVE CONTROL	(ppb)	SENSITIVITIES	(ppb)
Quantitative Assay [International Market]	30 minutes	Aflatoxin M1 [Dilute 1:9 to make 25 ppt standard] □	0.25	Aflatoxin M1	.025. 25
Semi-Quantitative [U.S. Market]	15 minutes	Aflatoxin M1	0.25	Aflatoxin M1	Confirm at 500 ppt

