


# PocketSwab Plus ATP Swab for Surface Hygiene Verification

**5  
seconds**

## Product Overview



The Charm® **PocketSwab® Plus ATP Swab (PocketSwab Plus swab)** is a room-temperature stable, self-contained, single service swab that enables total surface hygiene verification for ATP (adenosine triphosphate) in only 5 seconds, allowing real-time corrective action.


This high sensitivity swab is relied on for cleaning validation, prevention of cross-contamination, infection control, HACCP audits and enforcement across multiple industries.

### Industries used for:



Animal Care	Beverage	Clean Room	Cosmetics	Environmental
Food	Grain	Health Care	Medical Devices	Nutraceuticals
Personal Care Products	Pharmaceuticals	Pulp & Paper	Truck Wash	Water Treatments
Wineries				

### Benefits of PocketSwab Plus Swab Include:


- ◆ **5 Second results = Faster Pre-op Inspection**
- ◆ **Enhanced sensitivity to ATP = Reduced microbial counts; Reduced consumer complaints**
- ◆ **Room temperature stable with 12-month shelf life = Reduced freight cost with bulk packaging**
- ◆ **Recyclable components = Reduced waste**



## Simple Procedure



## Swab Extender



More Charm customers are incorporating Swab Extenders to monitor hard-to-reach and inaccessible areas in their Sanitation Standard Operating Procedures (SSOPs).

Available in 3 and 16 foot (0.9 m and 4.9 m) lengths, users can simply lock the relevant Charm ATP swab into the extender and adjust the angle to reach their target.

## Ordering Info

Order Codes	Each Kit Includes
<b>POCK-PLUS-25K</b>	PocketSwab Plus 25 Kit
<b>POCK-PLUS-100K</b>	PocketSwab Plus 100 Kit
<b>POCK-PLUS-600K</b>	PocketSwab Plus 600 Kit
<b>POCK-PLUS-1000K</b>	PocketSwab Plus 1000 Kit
<b>150-0035</b>	3 ft (0.9 m) Extender
<b>150-0040</b>	16 ft (4.9 m) Extender

See Operator's Manual for additional kits and order codes