



ATP Monitoring Program for the Hospitality Industry

Sanitation Verification Simplified

Charm's trusted adenosine triphosphate (ATP) detection technology and sophisticated novaLINK software is a powerful system that collects data to evaluate and improve your cleanliness program.

The novaLUM® II-X system wirelessly syncs data to local or corporate databases. Automated reports can be emailed to Operation Managers for sanitation verification.



novaLUM II-X system



PocketSwab® Plus ATP Swab

What is Adenosine Triphosphate (ATP)?

An energy-carrying molecule found in living or once-living cells.

Can the ATP tests measure for viruses?

While it doesn't measure for viruses, ATP surface testing specifically identifies the presence or absence of biological matter, which is a gauge of general sanitation and levels of hygiene.



Quick and Easy Surface Testing



Certified Inspectors



Sync Results and Reports



Inspected Clean

www.charm.com

T +1.978.687.9200 | F +1.978.687.9216 | E info@charm.com | 659 Andover Street, Lawrence, MA 01843-1032, USA
© 2021 Charm Sciences, Inc. | Charm, novaLUM, and PocketSwab are registered trademarks of Charm Sciences, Inc. | MRK-1434-001 Jul-2021

