



## Mycotoxin Charm® Case

The mycotoxin charm case is a lab-on-the-go with everything you need to make mycotoxin testing easy



**The Mycotoxin Charm Case is customizable with the following options:**

- Charm EZ®-M System
- ROSA® Quad Incubator
- Printer (optional)
- Centrifuge and microtubes
- Pipettes
- Whirlpack bags
- Graduated cylinder
- Canister of ROSA test strips

**The Charm Myco Case fits all the necessary equipment and consumables into a convenient carrying case.**

- Designed for maximum portability – easy to transport and setup
- Fully customizable
- Can be used with any of Charm's ROSA mycotoxin tests

**Scan the QR code to speak to our mycotoxin experts!**

