

INTRODUCING

Peel Plate® INTELLICOUNT System

Automated colony counting for
Peel Plate microbial tests

With the Peel Plate IntelliCount System:



Cut training time with an intuitive, automated workflow



Increase throughput with fast, standardized counts



Minimize errors with patented dome lighting and SphereFlash analysis




Maintain audit-ready digital records (Excel, PDF, LIMS export)



Choose the right software: PRO for advanced compliance, AI for operator-free accuracy (computer with internet connection required)




Fast, accurate, and consistent colony counting designed for busy microbiology labs. The IntelliCount System automates colony counting with adapters for both 1 mL and 5 mL Peel Plate tests. The system connects to your computer, running SphereFlash® software for instant, traceable results.




Count

Place your developed Peel Plate in the holder, press start, and the dome lowers to provide precise lighting without ambient light influences.



Capture

SphereFlash software processes colonies and stores image results and sample information.



Refine

Operators can fine-tune with adjustable settings (contrast, size, colony spread). With SphereFlash software options, you choose the level of automation.

www.charm.com

T +1.978.687.9200 | F +1.978.687.9216 | E info@charm.com | 659 Andover Street, Lawrence, MA 01843-1032

© 2025 Charm Sciences, Inc. | Charm, IntelliCount, and Peel Plate are registered trademarks of Charm Sciences, Inc. | MRK-1606-001 Sep-2025

CHARM
SCIENCES INC.
NOTHING WORKS LIKE A CHARM